



Product Overview

PRODUCT:	Alpha-1 Acid Glycoprotein Elevated, Defibrinated Human Plasma
PART NUMBER:	4819
FORM / APPEARANCE:	Liquid, amber to straw, free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	No preservatives
STORAGE:	≤ - 20°C
EXPIRATION:	2 years from DOM

Analytical Properties

ANALYTE	VALUE
TOTAL PROTEIN:	7 – 11 g/dL
ALPHA-1 ACID GLYCOPROTEIN:	> 250 mg/dL
TURBIDITY (O.D. @ 546nm):	< 0.5
pH:	7.0 – 7.5
BIOBURDEN:	< 100 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS