



## Product Overview

PRODUCT:	Bone and Anemia Control, 2 Level
PART NUMBER:	4905
FORM / APPEARANCE:	Off-white crystalline powder
FILTRATION:	0.22 micron-filtered
ADDITIVES:	No preservatives
STORAGE:	2 – 8 °C
EXPIRATION:	3 years from date of manufacture

## Analytical Properties

ANALYTE	VALUES	
	Level 1	Level 2
PTH:	55 – 65 pg/mL	1425 – 1575 pg/mL
VITAMIN D:	15 – 25 ng/mL	65 – 75 ng/mL
CALCITONIN:	95 – 105 pg/mL	595 – 605 pg/mL
FOLATE:	2 – 6 ng/mL	12 – 18 ng/mL
VITAMIN B12:	245 – 255 pg/mL	950 – 1050 pg/mL
BIOBURDEN:	< 100 CFUs/mL	< 100 CFUs/mL

### **FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY**

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS