



Product Overview

PRODUCT:	Bovine Cholesterol Concentrate
PART NUMBER:	5334
FORM / APPEARANCE:	Liquid, amber, free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	None
STORAGE:	≤ -40 °C
EXPIRATION:	2 years from date of manufacture at ≤ -40°C 1 year from date of manufacture at ≤ -20°C

Analytical Properties

ANALYTE	VALUE
TOTAL CHOLESTEROL:	≥ 3000 mg/dL
TOTAL PROTEIN:	< 8.0 g/dL
TRIGLYCERIDES:	< 500 mg/dL
CALCIUM:	< 5 mg/dL
SODIUM:	< 70 mmol/L
CHLORIDE:	< 60 mmol/L
POTASSIUM:	< 5 mmol/L
AMMONIA:	< 1.0 mg/dL
pH:	7.8 – 8.2
BIOBURDEN:	< 100 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

Prepared from bovine serum obtained from cattle less than 30 months of age and that have received ante and post mortem health inspections from US Veterinary Service Inspectors within facilities inspected by the United States Department of Agriculture.