



Product Overview

PRODUCT:	Bovine Serum Albumin
PART NUMBER:	5298
FORM / APPEARANCE:	Off-white crystalline powder
FILTRATION:	0.22 micron-filtered
ADDITIVES:	None
STORAGE:	2 – 8 °C
EXPIRATION:	3 years from date of manufacture, reassign after 3 years

Analytical Properties

ANALYTE	VALUE
TOTAL PROTEIN:	> 96% by Biuret
ALBUMIN PURITY:	> 98% by Electrophoresis
BIOBURDEN:	< 10 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

Prepared from bovine serum obtained from cattle less than 30 months of age and that have received ante and post mortem health inspections from US Veterinary Service Inspectors within facilities inspected by the United States Department of Agriculture.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.