



Product Overview

PRODUCT:	Fibrinogen Deficient Plasma, Defibrinated
PART NUMBER:	5422
FORM / APPEARANCE:	Liquid, amber to straw, free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	20 mM HEPES
STORAGE:	≤ - 60 °C
EXPIRATION:	1 year from date of manufacture; Stable when thawed for 4 hours at 2 – 8 °C

Analytical Properties

ANALYTE	VALUE
FIBRINOGEN:	1 g/dL
pH:	7.0 – 7.8
BIOBURDEN:	< 100 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS