



## Product Overview

PRODUCT:	Human Serum Albumin, 25% Solution
PART NUMBER:	3447
FORM / APPEARANCE:	Liquid, clear to straw, free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	No preservatives
STORAGE:	2 – 8 °C
EXPIRATION:	5 years from date of manufacture

## Analytical Properties

ANALYTE	VALUE
TOTAL PROTEIN:	24 – 26 g/dL
ALBUMIN PURITY:	> 96 %
CALCIUM:	< 1 mg/dL
SODIUM:	< 200 mEq/L
CHLORIDE:	< 100 mEq/L
POTASSIUM:	< 2 mEq/L
pH:	6.7 – 7.3
BIOBURDEN:	< 10 CFUs/mL

### **FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY**

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according of FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS