



Product Overview

PRODUCT:	Off the Clot, Normal Human Serum
PART NUMBER:	2760
FORM / APPEARANCE:	Liquid, clear yellow, free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	None
STORAGE:	<-20 °C
EXPIRATION:	5 years from date of manufacture

Analytical Properties

ANALYTE	VALUE
TOTAL PROTEIN:	6.3-7.0
CALCIUM:	8.5-11.5
D-DIMER:	<1.2
TROPONIN I:	<0.50
FIBRINOGEN:	Negative
SODIUM:	130-148
POTASSIUM:	3.2-5.4
BICARBONATE:	10-30
CHLORIDE:	90-110
GLUCOSE:	50-120
BUN:	5-20
CREATININE:	0.4-1.3
TOTAL BILIRUBIN:	<1
ALBUMIN:	3.5-5.5
pH:	7.1-7.9
BIOBURDEN:	<10

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human blood units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.