



Product Overview

PRODUCT:	Total Human Cholesterol Concentrate
PART NUMBER:	4567
FORM / APPEARANCE:	Liquid, clear amber, free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	No preservatives
STORAGE:	< - 60 °C
EXPIRATION:	2 years from date of manufacture

Analytical Properties

ANALYTE	VALUE
TOTAL CHOLESTEROL:	> 2,500 mg/dL
TRIGLYCERIDE:	As Found
TRIG. RATIO (TRIG:TC):	< 0.8
TOTAL PROTEIN:	< 16 g/dL
pH:	7.5 – 8.5
BIOBURDEN:	< 10 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS