



Product Overview

PRODUCT:	Whole Blood Glucose Control, 2 Level
PART NUMBER:	5630
FORM / APPEARANCE:	Liquid
FILTRATION:	0.22 micron-filtered
ADDITIVES:	Stabilizers
STORAGE:	2 – 8 °C
EXPIRATION:	6 months from date of manufacture

Analytical Properties

ANALYTE	VALUES	
	Level 1	Level 2
GLUCOSE:	63 – 77 mg/dL	126 – 154 mg/dL
HEMOGLOBIN:	12 – 16 g/dL	12 – 16 g/dL
pH:	7.2 – 7.6	7.2 – 7.6
BIOBURDEN:	< 100 CFUs/mL	< 100 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

Individual donor units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according of FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.