



Product Overview

PRODUCT:	ASO Concentrate
PART NUMBER:	4333
FORM / APPEARANCE:	Liquid, Light brown to Amber, Free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	None
STORAGE:	≤-20°C
EXPIRATION:	3 years from date of manufacture

Analytical Properties

ANALYTE	VALUE
ASO:	4200 IU/mL
TOTAL PROTEIN:	As Found
IgG:	< 2000 mg/dL
ASO/T. PROTEIN RATIO:	>40 IU/mg
ASO/IgG RATIO:	>500 IU/mg
pH:	7.4
TURBIDITY:	< 0.5
BIOBURDEN:	< 100 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

Individual donor units in the preparation of this product have been tested and found to be non-reactive for HBsAg and negative for anti-HIV I/II, anti-HCV.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.