



## Product Overview

PRODUCT:	Ceruloplasmin, Partially Purified
PART NUMBER:	5010
FORM / APPEARANCE:	Liquid, blue, free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	No preservatives
STORAGE:	≤ - 20°C
EXPIRATION:	3 years from date of manufacture

## Analytical Properties

ANALYTE	VALUE
CERULOPLASMIN:	> 500 mg/dL
pH:	7.0 – 7.5
BIOBURDEN:	< 100 CFUs/mL

\*Testing performed at a 5% solution with sterile DI H<sub>2</sub>O

### **FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY**

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS