



## Product Overview

PRODUCT:	Human Platelet Lysate
PART NUMBER:	5015
FORM / APPEARANCE:	Liquid, amber to straw, product may appear cloudy and contain particulates after thaw
ADDITIVES:	None
STORAGE:	≤ - 20 °C
EXPIRATION:	5 years from date of manufacture

## Analytical Properties

ANALYTE	VALUE
MYCOPLASMA:	Negative
ENDOTOXIN:	< 0.5 EU/mL
STERILITY TESTING:	No Growth

### **FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY**

Individual donor units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human platelet units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.