



Product Overview

PRODUCT:	COVID-19 Antibody Negative Defibrinated Human Plasma from US Origin Source Plasma
PART NUMBER:	6095
FORM / APPEARANCE:	Liquid, straw to amber, free of visible foreign matter or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	None
STORAGE:	≤-10 °C
EXPIRATION:	5 years from date of manufacture

Analytical Properties

ANALYTE	VALUE
COVID-19 IgM Ab:	Negative
COVID-19 IgG Ab:	Negative
TOTAL PROTEIN:	6.0 – 6.5 g/dL
TOTAL CHOLESTEROL:	As Found
TRIGLYCERIDE:	As Found
CALCIUM:	As Found
SODIUM:	As Found
CHLORIDE:	As Found
POTASSIUM:	As Found
AMMONIA:	As Found
pH:	6 – 8
BIOBURDEN:	< 10 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human blood units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HBsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, HBV DNA and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.

