



Product Overview

PRODUCT:	COVID-19 Antibody Total Negative (IgG+IgM) Human Plasma
PART NUMBER:	6101
FORM / APPEARANCE:	Liquid, straw to amber, free of visible foreign matter or precipitates
FILTRATION:	Available
ADDITIVES:	None
STORAGE:	≤-20 °C
EXPIRATION:	5 years from date of manufacture

Analytical Properties

ANALYTE	VALUE
COVID-19 IgM Ab:	Negative
COVID-19 IgG Ab:	Negative
TOTAL PROTEIN:	As Found
SODIUM:	As Found
POTASSIUM:	As Found
CHLORIDE:	As Found
pH:	7.1-7.5
BIOBURDEN:	< 100 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human blood units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HBsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, HBV DNA and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.