



Product Overview

PRODUCT:	Heat Inactivated Donor Goat Serum
PART NUMBER:	6100
FORM / APPEARANCE:	Liquid, straw to amber, clear and odorless. May have some precipitates upon thawing
FILTRATION:	0.22 micron filtered
ADDITIVES:	1.0 g/L Sodium Azide
STORAGE:	≤ -20 °C
EXPIRATION:	5 years from date of manufacture

Analytical Properties

ANALYTE	VALUE
TOTAL PROTEIN:	60 – 90 mg/mL
IgG:	As Found
OSMOLALITY:	As Found
pH:	As Found
HEMOGLOBIN:	< 30 mg/dL
STERILITY:	Negative for Growth

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

Obtained from healthy donor goats that have received clean health inspections from licensed veterinarians in good standing in their respective country, and the animals are of specified country's origin only.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.