



## Product Overview

PRODUCT:	Human Hemoglobin Solution, < 5%
PART NUMBER:	4635
FORM / APPEARANCE:	Liquid
FILTRATION:	None
ADDITIVES:	Stabilizers added
STORAGE:	≤-20°C , protected from light
EXPIRATION:	2 years from date of manufacture

## Analytical Properties

ANALYTE	VALUE
HEMOGLOBIN:	< 5%
PH:	6.5 – 7.5

### **FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY**

All human blood units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.