



Product Overview

PRODUCT:	ASO/RF/CRP Control, 2 Level
PART NUMBER:	5324
FORM / APPEARANCE:	Liquid, amber to straw, free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
STORAGE:	≤ - 20 °C
EXPIRATION:	2 years from date of manufacture

Analytical Properties

ANALYTE	VALUES	
	Level 1	Level 2
ASO:	500 – 600 IU/mL	500 – 600 IU/mL
RF:	150 – 200 IU/mL	150 – 200 IU/mL
CRP:	20 – 30 mg/L	80 – 100 mg/L
TOTAL PROTEIN:	5 – 8 g/dL	5 – 8 g/dL
TURBIDITY (O.D. @ 546nm):	< 0.8	< 0.8
BIOBURDEN:	< 100 CFUs/mL	< 100 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS