



Product Overview

PRODUCT:	C4 Concentrate
PART NUMBER:	5321
FORM / APPEARANCE:	Liquid, amber to straw, free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	No preservatives
STORAGE:	<-20 °C
EXPIRATION:	2 years from date of manufacture

Analytical Properties

ANALYTE	VALUE
TOTAL PROTEIN:	6-11 g/dL
C4 CONCENTRATE:	> 70
pH:	7.6 – 8.1
TURBIDITY:O.D. @546nm	As Found
BIOBURDEN:	< 100 CFUs/mL

*Testing performed at a 5% solution with sterile DI H₂O

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS