



Product Overview

PRODUCT:	Elevated A1c Red Blood Cells
PART NUMBER:	3223
FORM / APPEARANCE:	Liquid, Dark, cherry red with no clots or abnormal particulate
FILTRATION:	Filtration available
ADDITIVES:	K3 EDTA Anticoagulant
STORAGE:	2 - 8 °C
EXPIRATION:	2 months from date of manufacture

Analytical Properties

ANALYTE	VALUE
HbA1c:	> 10%

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human blood units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.