



## Product Overview

PRODUCT:	hCG Negative Urine, Unprocessed
PART NUMBER:	5439
FORM / APPEARANCE:	Liquid, light yellow
FILTRATION:	None
ADDITIVES:	None
STORAGE:	< - 20 °C
EXPIRATION:	3 years from date of manufacture

## Analytical Properties

ANALYTE	VALUE
MICROALBUMIN:	< 0.5 mg/dL
SPECIFIC GRAVITY:	≥ 1.010
GLUCOSE:	Negative
NITRITE:	Negative
PROTEIN:	Negative
KETONE:	Negative
BILIRUBIN:	Negative
UROBILINOGEN:	0.2 – 1.0 EU/dL
LEUKOCYTE:	Negative
BLOOD:	Negative
hCG:	Negative
pH:	As Found

### **FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY**

Donors have been tested by FDA approved tests and found to be negative or non-reactive for HBsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, HCV RNA. Donors were tested and found negative for RPR according to FDA guidelines.

**CAUTION: FOR LABORATORY REAGENT USE OR MANUFACTURING INTO NONINJECTABLE PRODUCTS ONLY. THE USE OF THIS MATERIAL IN HUMAN MEDICATION IS PROHIBITED.**