



Product Overview

PRODUCT:	anti-HBs, anti-HBc, anti-TP Negative Off Clot Human Serum
PART NUMBER:	5606
FORM / APPEARANCE:	Liquid, amber to straw, free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	None
STORAGE:	≤ - 10 °C
EXPIRATION:	5 years from date of manufacture

Analytical Properties

ANALYTE	VALUE
TOTAL PROTEIN:	≥ 6.0 g/dL
aHBs:	< 0.5 mIU/mL
aHBc:	< 0.3 S/CO
aTP:	Negative
AMMONIA:	< 0.5 mg/dL
pH:	7.0 – 7.5
TURBIDITY (O.D. @ 546nm):	< 0.16
BIOBURDEN:	< 100 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according of FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.