



## Product Overview

PRODUCT:	CRP Control, Liquid, 2 Level
PART NUMBER:	5577
FORM / APPEARANCE:	Liquid, amber to straw, free of visible foreign matter, clots, or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	Stabilizers
STORAGE:	2 – 8 °C
EXPIRATION:	2 years from date of manufacture

## Analytical Properties

ANALYTE	VALUES	
	Level 1	Level 2
CRP:	17 – 23 mg/L	90 – 110 mg/L
TOTAL PROTEIN:	5.0 – 6.0 g/dL	5.0 – 6.0 g/dL
pH:	7.2 – 7.6	7.2 – 7.6
BIOBURDEN:	< 100 CFUs/mL	< 100 CFUs/mL

### **FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY**

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

**CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.**

**CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS**