



Product Overview

PRODUCT:	Lipid Linearity, Liquid, 6 Levels
PART NUMBER:	5648
FORM / APPEARANCE:	Liquid, Clear amber, free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	Proclin 300
STORAGE:	≤ -20°C
EXPIRATION:	1 year from date of manufacture

Analytical Properties

ANALYTE	VALUES					
	Level 1	Level 2	Level 3	Level 4	Level 5	Level 6
TOTAL CHOLESTEROL (mg/dL):	90 – 110	174 – 212	257 – 315	342 – 418	426 – 520	509 – 623
TRIGLYCERIDES (mg/dL):	90 – 110	168 – 206	246 – 300	324 – 396	401 – 491	480 – 586
HDL CHOLESTEROL (mg/dL):	31 – 39	51 – 63	71 – 87	91 – 111	111 – 135	131 – 160
LDL CHOLESTEROL (mg/dL):	As Found	As Found	As Found	As Found	As Found	As Found
TOTAL PROTEIN (g/dL):	4.5 – 10	4.5 – 10	4.5 – 10	4.5 – 10	4.5 – 10	4.5 – 10
pH:	7.6 – 8.0	7.6 – 8.0	7.6 – 8.0	7.6 – 8.0	7.6 – 8.0	7.6 – 8.0
BIOBURDEN (CFUs/mL):	< 100	< 100	< 100	< 100	< 100	< 100

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.