



Product Overview

PRODUCT:	Porcine LDL Cholesterol Concentrate
PART NUMBER:	5762
SOURCE:	Porcine Serum certified as U.S. Origin
FORM / APPEARANCE:	Liquid, Opaque, Yellow to straw, free of visible foreign matter, clots, or precipitate
ELECTROPHORESIS:	One major band corresponding to LDL
FILTRATION:	0.22 micron-filtered
STORAGE:	2 – 8 °C
EXPIRATION:	1 year from date of manufacture

Analytical Properties

ANALYTE	VALUE
TOTAL CHOLESTEROL:	> 3000 mg/dL
TOTAL PROTEIN:	< 8.0 g/dL
IRON:	< 200 µg/dL
CALCIUM:	< 5 mg/dL
SODIUM:	< 150 mEq/L
CHLORIDE:	< 150 mEq/L
AMMONIA:	< 400 µmol/L
BIOBURDEN:	< 10 CFUs/mL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

Prepared from porcine serum collected from healthy swine at licensed facilities inspected by the United States Department of Agriculture. Collection was performed in accordance with applicable regulations.