



## Product Overview

PRODUCT:	ASO Elevated, Delipidized, Defibrinated Human Plasma
PART NUMBER:	4866
FORM / APPEARANCE:	Liquid, Straw to Amber, Free of visible foreign matter, clots or precipitates
FILTRATION:	0.22 micron-filtered
ADDITIVES:	None
STORAGE:	≤-20°C
EXPIRATION:	5 years from date of manufacture

## Analytical Properties

ANALYTE	VALUE
ASO:	>1200 IU/mL
TOTAL PROTEIN:	< 9 g/dL
TRIGLYCERIDES:	< 20 mg/dL
CHOLESTEROL:	< 20 mg/dL
pH:	6.8 – 7.2
TURBIDITY:	< 0.8
BIOBURDEN:	< 100 CFUs/mL

### **FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY**

Individual donor units in the preparation of this product have been tested and found to be non-reactive for HBsAg and negative for Anti-HIV I/II, anti-HCV, HIV-I RNA and HCV RNA. Donors of human donor units used in the making of this product were tested and found negative for Syphilis.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.